



Spherical Lesion Phantom

GAMMEX4 08 LE

Ultrasound image quality in three dimensions: Axial, lateral and elevational.

The Spherical Lesion Phantom Gammex4 08 LE provides a unique way of testing resolution performance of ultrasound scanners. The Spherical Lesion Phantom contains 2 mm and 4 mm diameter tissue mimicking spherical lesions which lie on a single plane at the center of the phantom. Axial, lateral and elevational resolution are accounted for simultaneously and equally for all types of ultrasound systems and configurations. In the 2 mm section, there are 105 anechoic spheres at 0.5 cm depth intervals and in the 4 mm section there are 211 anechoic spheres at 0.75 cm depth intervals.

The Spherical Lesion Phantom Gammex4 08 LE incorporates our Tissue Mimicking (TM) gel which provides a

smoother background texture than conventional tissue mimicking gels. The Gammex4 08 LE allows production of lesions with negligible echogenicity while producing no distal enhancement or shadowing inherent with other gels. The TM gel is also optimized for use with tissue harmonics imaging technology.

In addition, the Gammex4 08 LE has an new composite film scanning surface that has improved transmission properties so more of the ultrasonic beam can be transmitted and received.

Optional accessories for the Gammex4 08 LE include Soft Foam-Lined Carrying Case, Rigid Case, or the Precision Ultrasound Transducer Guide.

continued

WWW.GAMMEX.COM

PEO B.V.

info@gotopeo.com
www.gotopeo.com

The Netherlands

Havenweg 16, 6603 AS Wijchen
+31 (0)24 648 86 88

Belgium

Watermolenstraat 2, B-2910 Essen
+32 (0)3 309 32 09

CoC 34107894

VAT NL807859151B01

IBAN NL29 RABO 0356 1960 46

BIC RABONL2U






GAMME X4 08 LE

ULTRASOUND

continued from front...



This ultrasound image demonstrates the anechoic spheres contained in the Spherical Lesion Phantom Gammex4 08 LE. Note the well defined 2 mm and 4 mm targets.

SPECIFICATIONS

High Resolution Tissue Mimicking Gel

Speed of sound .. 1540 ± 10 m/s at 22 °C

Attenuation 0.50 ± 0.07 ± 0.05 dB/cm/MHz

Anechoic Spherical Lesions

Size..... 4 mm in ap plane 0.5 to 6 cm

2 mm in a plane 0.5 to 10.5 cm

Contrast..... -30 dB relative to background

Construction

Scanning

surface Composite film

Walls Extruded ABS

Dimensions 23.2 x 8.25 x 18.5 cm
(9.25 x 3.25 x 7.25 in)

Weight..... 2.8 kg (6 lbs 5 oz)

Acoustic measurements at 4.5 MHz, 22 °C

PEO B.V.

info@gotopeo.com
www.gotopeo.com

The Netherlands

Havenweg 16, 6603 AS Wijchen
+31 (0)24 648 86 88

Belgium

Watermolenstraat 2, B-2910 Essen
+32 (0)3 309 32 09

CoC 34107894

VAT NL807859151B01

IBAN NL29 RABO 0356 1960 46

BIC RABONL2U

