



DR/CR/Fluoroscopy Solutions

Density Control Function (DCF) Test Tool

PN 804964



Base Plate:	Aluminum
Exposure Plate	Stainless Steel
Overall	15.2 x 30.5 x 0.95 cm (HWD)
Window	12 x 18.6 mm (0.47 x 1.12 in)
Weight	0.8 kg (1.75 lbs)
File Size	7 x 9.5 in
Exposure Steps	11 (-1 to -5, zero, +1 to +5)
Compatibility	Tissue Equivalent Breast material
Notes	Tissue Equivalent Material sold separately

PEO B.V.

info@gotopeo.com
www.gotopeo.com

The Netherlands

Havenweg 16, 6603 AS Wijchen
+31 (0)24 648 86 88

Belgium

Watermolenstraat 2, B-2910 Essen
+32 (0)3 309 32 09

CoC 34107894

VAT NL807859151B01

IBAN NL29 RABO 0356 1960 46

BIC RABONL2U

